# TERUMO® SurGuard®3 CE 0197

**Non-toxic. Non-pyrogenic.**

**The device has no components made of natural rubber latex.**

## SAFETY HYPODERMIC NEEDLE

The TERUMO® SurGuard®3 SAFETY HYPODERMIC NEEDLE is intended for use in the aspiration and injection of fluids for medical purposes. The Terumo Safety Hypodermic Needle is compatible for use with standard luer slip and luer lock syringes. Additionally, after withdrawal of the needle from the body, the attached needle safety sheath can be manually activated to cover the needle immediately after use to minimize risk of accidental needlestick.


## HYPODERMIC SYRINGE WITH SAFETY NEEDLE

The TERUMO® SurGuard®3 HYPODERMIC SYRINGE WITH SAFETY NEEDLE is intended for use in the aspiration and injection of fluids for medical purposes. Additionally, after withdrawal of the needle from the body, the attached needle safety sheath can be manually activated to cover the needle immediately after use to minimize risk of accidental needlestick.

## INSTRUCTIONS FOR USE

### STEP 1.


Tighten the syringe to the needle using aseptic technique. **Grip the base of the needle**, not the safety sheath, push and turn the syringe clockwise.

### STEP 2.

Move the safety sheath away from the needle and toward the syringe barrel **to the angle shown**, prior to removing the needle cap. Perform injection procedure according to established technique.


(Note: The needle "bevel up" position is oriented to the position of the safety sheath.)

### STEP 3.


After the injection procedure, **use a one-handed technique to activate the safety mechanism** using any of the three (3) methods illustrated below.

(Activation is verified by an audible and/or tactile "click", and can be visually confirmed.)

### STEP 4.

Visually confirm that the safety sheath is fully engaged.

### STEP 5.

Dispose of used needles and materials following the policies and procedures of your facility as well as federal and local regulations for sharps disposal.

## WARNINGS

- HANDLE WITH CARE TO AVOID NEEDLESTICKS.

Dispose of used needles and materials following the policies and procedures of your facility as well as federal and local regulations for sharps disposal.

## WARNINGS

- HANDLE WITH CARE TO AVOID NEEDLESTICKS.
- DO NOT REPLACE ORIGINAL NEEDLE COVER AFTER USE.
- USE ONCE AND DISCARD IMMEDIATELY IN ACCORDANCE WITH LOCAL SAFETY STANDARDS.

## CAUTIONS

- If needle is bent or damaged, no attempt should be made to straighten needle or use the product.
- Do not attempt to deactivate the safety device by forcing the needle out of the safety sheath.
- For single use only. Do not reuse. Do not resterilize.

## PRECAUTIONS

- Keep hands behind needle at all times during use and disposal.
- Observe Universal Precautions on ALL patients.
- Do not use if the unit package or product has been damaged or contaminated.
- Do not store at extreme temperature and humidity. Avoid direct sunlight.

 Manufacturer TERUMO (PHILIPPINES) CORPORATION  
124 EAST MAIN AVE., LAGUNA TECHNOPARK, BINAN, LAGUNA, PHILIPPINES  
 TERUMO EUROPE N.V., Interleuvenlaan 40, 3001 LEUVEN, BELGIUM  
TERUMO MEDICAL CORPORATION, 2101 Cottontail Lane, SOMERSET, NJ 08873 USA

	Do not reuse
	Consult instructions for use
<b>Rx ONLY</b>	CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Sterilized using irradiation (electron beam)
	Batch code
	Use by date
	Catalogue number
	Manufacturer
	Authorized representative in the european community
	Do not use if package is damaged

 and TERUMO are registered trademarks of Terumo Corporation, Tokyo, Japan and they are registered with the U.S. Patent & Trademark Office.

**SurGuard<sup>®3</sup>** is a registered trademark of Terumo Medical Corporation.  
01/2013 REV.0

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